

# **EASYBIL-P®**

# Bilirubin Analyzer

Advanced and Optimized instrument for measuring Neonate's Bilirubin









#### **Main Features**

- EasyBil-P provides quick, easy and accurate results of total bilirubin in whole blood of pediatric samples with volume of 50-80  $\mu$ L
- Eliminates interference from turbidity and hemolysis with its unique dual wavelength measurement system
- · In-built self diagnostics check
- Life long LED lighting source that does not require replacement
- New design with auto locking capillary holder that prevents damage & breakage of capillary tubes
- Quick reading photometer display results in 2 seconds
- Capillary based system, no reagents required



## **Specifications**

Light source: LED's

Wave Length: 460nm & 550nm

Measuring Range: 0-40 mg/dL

Display: Large display

Accuracy: Within 3.5% for values b/w 3.0 - 40 mg/dL

Within 5% for values b/w 0.0 - 3.0 mg/dL

Power Supply: 220V - 50Hz

Dimention: 222mm X 180mm X 100mm

Printer: 56mm Thermal graphic 384 dots/line

### **Accessories**

External Adapter (12V DC): 1pc

Capillary Tube, Heparinized: 1 Box (100pcs)

Sealing Wax: 1pc
Capillary Holder: 1pcs
User's Manual: 1pc
Fuses: 2pcs
Printer roll: 2pcs

Technical specifications subject to change without prior notice

Product manufactured by Microlab Instruments Pvt. Ltd.

CIN: U52609GJ2021PTC128192 An ISO 9001 & ISO 13485 Certified Company



Factory: 246/P1, K.T. Compound, Uma Estate, Nr. Bhagyoday Hotel,