



INCUBATOR

INCUBATOR series

Dry Bath Incubator



INCUBATOR - 30S



INCUBATOR - 30



INCUBATOR - 30D - Dual block



Incubator

Default Temperature: 37°C (±0.2°C)
Default Time: 3 min / 5 min /
10 min / 15min
Temperature Setting: 25°C to 60°C
3 Inbuilt Timers: T1 / T2 / T3
1 Location Open: Editable up to
99 hours / 59 min

Power / Dimensions

Voltage: 110 - 220V / 50 - 60 Hz
Power Supply: <60W
Temperature: 10°C to 45°C
Relative humidity: ≤80%
Dimensions: 22cm X 27cm X 21cm

Interface

Screen: Alpha Numeric 20X2 lines
Keyboard: Tactile Soft - Touch

Technical specifications subject to change without prior notice

Product manufactured by Microlab Instruments Pvt. Ltd.

CIN : U52609GJ2021PTC128192
An ISO 9001 & ISO 13485 Certified Company



Microlab Instruments Private Limited

Head Office: 3, Akshay Apt., Nr. "SANSKRUT", B/H. Popular House, Ashram Road, Ahmedabad - 380 009. India.

T: +91 79-2657 4862, +91 99798 79096

Factory: 246/P1, K.T. Compound, Uma Estate, Nr. Bhagyoday Hotel, Nr. Toll Plaza, Vasna-Iyava, Sanand, Ahmedabad - 382 170. India.

T: +91 079-2990 0453

E: info@microlab.co.in | W: www.microlab.co.in